



# Suction Bag and Accessories

- ✓ SAFE
- ✓ SIMPLE
- ✓ USE FRIENDLY



LILY MEDICAL

GMP, ISO 13485, ISO 14001, ISO 45001, ISO 17025 Certifications

## FEATURES



- 1 Serial Tube
- 2 Pour Spout
- 3 Patient Port
- 4 Handle
- 5 Vacuum Port



- 1 Unidirectional Valve



- 1 Pour Spout
- 2 Patient Port
- 3 Handle



- 1 Hydrophobic Suction Filter
  - >> Overflow protection
  - >> High air flow rate
  - >> Super Bacterial Filtration Efficiency (ASTM F2101-07) BFE 99.9%
- 2 Unidirectional Valve
- 3 VACUUM Port  
Inline with negative pressure system to provide suction in the bag

## SPECIFICATIONS

Model No.	Description
<b>Suction Bag</b>	
CSB-G105	1.0L, w/ Hydrophobic filter, 2 Holes
CSB-G155	1.5L, w/ Hydrophobic filter, 2 Holes
CSB-G205	2.0L, w/ Hydrophobic filter, 2 Holes
CSB-G10C5	1.0L, w/ Hydrophobic filter, 3 Holes
CSB-G15C5	1.5L, w/ Hydrophobic filter, 3 Holes
CSB-G20C5	2.0L, w/ Hydrophobic filter, 3 Holes
CSB-W205	2.0L, w/o Hydrophobic filter, 2 Holes
CSB-W20C5	2.0L, w/o Hydrophobic filter, 3 Holes
<b>Canister</b>	
CSB-C101	1.0L, Green
CSB-C151	1.5L, Green
CSB-C201	2.0L, Green or White
<b>Trolley</b>	
CSB-C20H12	2L, 1 gauge, 2-canister, H:75cm
CSB-C20H14	2L, 1 gauge, 4-canister, H:75cm
CSB-C20H24	2L, 2 gauge, 4-canister, H:75cm
CSB-C20L12	2L, 1 gauge, 2-canister, H:55cm
CSB-C20L14	2L, 1 gauge, 4-canister, H:55cm
CSB-C20L24	2L, 2 gauge, 4-canister, H:55cm
<b>Connection Tube</b>	
CSC-XXX	45cm, 180cm, 200cm, 230cm, 360cm Customized length



Suction Bag



Canister



Trolley

## ACCESSORIES



CSB-A01  
Regulator Module



CSB-A02  
Suction Gauge,  
0~700mmHg



CSB-A03  
Shut Off Valve




CSB-A04  
Wall Mount,  
Green or White




Connection  
Tube

## Suction liner linkage

### Single liner



### Serial Connection



#### Flexible

Bag and canister can be used individually or in a serial connection up to 4 bottles

#### Safe

Completely sealed for endurance

#### User-friendly

Error-free, using color (white or green) for easy identification

#### Canister

Graduation lines are easy to read and minimum graduation is 50 ml



GMP Certificate



ISO 13485:2016



### LILY MEDICAL CORPORATION

TEL +886-37-612-625  
 FAX +886-37-612-720  
 E-MAIL [customer.service@lily-medical.com](mailto:customer.service@lily-medical.com)  
 WEB [www.lily-medical.com](http://www.lily-medical.com)  
 ADD No.28-2, Shun Jeau Diann, Chunan Town, Miaoli County, 35056, Taiwan