



LILY Silicone Nasogastric Tube

✔ Smoothness ✔ Comfort ✔ Longevity

Enhanced insertion smoothness, superior to PVC materials.
Comfortable and hypoallergenic, reducing patient discomfort.



LILY MEDICAL

Features

- + Made of 100% medical grade silicone
- + Available in one-piece or two-piece design
- + Two-piece design allows removal of outer tube when not feeding
- + Straight and smooth connection between inner and outer tubes prevents nutritional residue buildup and bacterial growth
- + Corrosion resistant to reduce frequency of replacement
- + Moderate flexibility for smooth insertion without damaging mucosa, reducing discomfort during long-term placement
- + X-ray opaque line
- + Graduated markings for easy depth measurement

Image



One-piece



Two-piece

Specification

Name	Model	Spec.
One-piece Silicone Nasogastric Tube	MDST-1-16FR120	16FR, Length 120cm
One-piece Silicone Nasogastric Tube	MDST-1-18FR120	18FR, Length 120cm
Two-piece Silicone Nasogastric Tube	MDST-2-16FR120	16FR, Length 60+60cm

This product is valid within three years

STERILE EO



Single use



LILY MEDICAL

TEL +886-37-612-625
FAX +886-37-612-720
E-MAIL customer.service@lily-medical.com
WEB www.lily-medical.com
ADD NO.28-2, Shun Jeau Diann, Chunan Town.
Miaoli County, 35056. Taiwan R.O.C

